JAJ SOLUTIONS

SKU : JAJ0053

EAN : 3701801600525

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User manual

Manual de uso

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Manuale d'uso

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DESCRIPTION

PRÉ-REQUIS

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MANUFACTURER

GUANGZHOU DONGDI TECHNOLOGY CO. LTD
202, BUILDING C, WU HAO, LONG QUAN XI LU,
SHI LING ZHEN, HUADU DISTRICT, GUANGZHOU CITY,
GUANGDONG PROVINCE, CHINA
MODEL NUMBER : AX188 / JAJ0053

DESCRIPTION

PREREQUISITES

Before first use, make sure to fully charge your mask in order to extend battery life.

DESCRIPTION

An LED facial mask uses different light colors to act on the skin without pain or harsh UV exposure. This is called LED light therapy.

The most common uses are:

Red light: stimulates collagen production and helps reduce wrinkles and improve firmness.

Blue light: helps reduce acne-causing bacteria.

Green light: may help reduce pigmentation spots and improve skin tone uniformity.

Yellow light: used to soothe redness and sensitive skin.

White light: used for skin regeneration and recovery.

USE

START-UP

- Adjust the mask to your face using the strap.
- Press the power button for 2-3 seconds.
- Select the desired color/light according to your needs.

DURING THE SESSION

- Keep your eyes closed.
- Relax for the recommended duration (15 minutes).
- Do not use for longer than recommended.

END OF THE SESSION

- Turn off the device.
- Remove the mask.
- Optionally apply a serum or moisturizer.
- Clean the mask with a soft cloth.

TERMS OF USE

BEFORE USE

- Read the entire manual before use.
- Use only on clean and dry skin.
- Remove makeup, thick cream, and jewelry.
- Check that the device is not damaged.
- Charge the device using only the recommended cable.

NORMAL USE

- Use only on the face (or indicated areas).
- Keep eyes closed during the session.
- Respect the recommended session duration (15 min).
- Follow the recommended frequency (e.g., 3-5 times per week).
- Do not use while sleeping.
- Use in a dry and well-ventilated environment.
- Stop using the device in case of discomfort, burning, or irritation.

TERMS OF USE

CONTRAINDICATIONS

Do not use:

- in case of photosensitive epilepsy,
- during pregnancy (according to manufacturer's recommendation),
- on damaged, burned, or infected skin,
- after a recent aggressive peel or laser treatment,
- in case of severe skin conditions without medical advice,
- if you are taking photosensitizing medications,
- on children without supervision.

EYE SAFETY

- Do not look directly at the LEDs.
- Use eye protection if provided.
- Stop immediately in case of visual discomfort.

ENVIRONMENTAL CONDITIONS

- Do not use in the rain or in a humid bathroom.
- Keep away from water and high temperatures.
- Do not expose to direct sunlight or fire.
- Recommended storage temperature: 0–40°C.

TERMS OF USE

BATTERY AND CHARGING

- Use only the recommended power supply.
- Do not disassemble, puncture, or heat the battery.
- Unplug after full charge.
- Do not use if the cable is damaged.

CLEANING AND MAINTENANCE

- Turn off before cleaning.
- Clean with a soft, slightly damp cloth.
- Do not immerse in water.
- Do not use harsh alcohol or solvents.

WARNINGS

- This product is intended for cosmetic use only.
- This product is not a medical device.
- Results may vary depending on the user.
- Keep out of reach of children.

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All trademarks are the property of their respective owners.
Importé par JAJ SOLUTIONS SAS, Fabriqué en Chine
Imported by JAJ SOLUTIONS SAS, Made in China.